SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCMSLQ

FINOFRAME COCR SOLID

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

	Gebrauchsanweisung	2
	Instructions for use	6
	Manuel d'utilisation	10
	Istruzioni per l'uso	14
	Instrucciones de uso	18
	Instrukcje użytkowania	22
	Instruktioner för användning	26
	Instruktioner til brug	30
	Bruksanvisning	34
	Käyttöohjeet	38
	Gebruiksaanwijzing	42
	Οδηγίες χρήσης	46
	Upute za uporabu	50
	Instrucțiuni de utilizare	54
	Naudojimo instrukcijos	58
	Návod k použití	62
	Návod na použitie	66
	Használati utasítás	70

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCMSLQ

FINOFRAME COCR SOLID

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

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Produkt/ Produkttyp

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Instructions for use

Product/ Product type

Sintered, hard Co/Cr blanks for the fabrication of individual dental restorations

Material type

Cobalt/chrome alloy

Class IIa medical device

Co/ Cr alloy type 4 according to DIN EN ISO 22674

The product is free from nickel, cadmium, beryllium and lead in accordance with DIN EN ISO 22674.

Users

Professional users who produce customised dental restorations

Target group

Patients of all ages and genders who require customised dental restorations

Indication / Intended use

The product is a medical device intended for the manufacture of dental prostheses for temporary or long-term use, which are partially introduced into the human body by means of a clinical procedure, but are inserted into the teeth or attached to the tooth structure with suitable fixing material and are therefore classified as Class IIa.

Indication:

- anatomically reduced crown and bridge frameworks in the anterior and posterior region
- Bridges with up to 14 units or bridges with small cross-sections
- Primary and secondary telescopic crowns
- clamps, bars, fasteners
- Fully anatomical crowns and bridges in the posterior region
- Restorations with thin cross-sections that are exposed to very high loads
- Free-end bridges with a maximum of one pontic
- implant-supported superstructures
- Removable partial dentures

The product is a type 4 Co/Cr alloy in accordance with DIN EN ISO 22674 and therefore has no indication restrictions compared to milled or cast type 4 Co/Cr alloys.

Contraindication

Do not use in case of proven hypersensitivity to one or more metals contained in the alloy

Veneering ceramics/ glaze

Co/Cr veneering ceramic

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCMSLQ

FINOFRAME COCR SOLID

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

Material properties/ Technical data

Composition:		Technical data (after final sintering):	
Co	59,5 - 69%	0.2% proof stress	475 MPa
Cr	26,5 - 30%	Percentage elongation at break	11%
Mo	4,5 - 7%	Percentage breakage constriction	14%
Mn	0 - 1%	E-modulus	186 GPa
Fe	0 - 1%	Tensile strength	731 MPa
Si	0 - 1%	Density	7.8 g/cm ³
C	0 - 0,35%	Corrosion resistance	< 200 µg/cm ²
other	< 1%	Tarnish resistance	Yes
		Vickers hardness	220 HV1
		Coefficient of thermal expansion	13.74 x 10 ⁻⁶ /K

Description

The product consists of blanks made of a biocompatible cobalt/chromium alloy for the direct fabrication of dental restorations. They are first axially pressed from an extremely fine powder and then individually isostatically compacted. They are then debinded and subsequently sintered to their final density and final hardness.

No further sintering process needs to be carried out by the user.

Processing

Scaffolding production

The desired restoration is milled from the selected blank. The blanks can be processed with all commercially available dental CAD/CAM milling machines.

Disconnect

Once the milling process is complete, the objects are cut out of the blank.

Cut out and grind frameworks and individual links using suitable carbide cutters or cutting discs. Always use the same rotating instruments for one alloy to avoid contamination.

Firing the ceramic

In principle, all commercially available ceramic bonding materials can be used. Please observe the processing instructions of the respective ceramic manufacturer and the CTE range specified therein for compatibility.

The minimum thickness of the finished copings should not be less than 0.3 mm. It is recommended to blast and clean (steam off) the frameworks with at least 110 µm aluminum oxide at 3-4 bar; the oxide firing is not mandatory, but can optionally be carried out for 5 minutes at 980 °C under vacuum (cleaning firing). The framework must be blasted with 110 µm aluminum oxide at 3-4 bar to thoroughly remove the oxide layer. Then clean the framework with a steam jet. When using a ceramic bonder, please follow the processing instructions of the respective manufacturer.

Soldering

We recommend a cobalt-based solder for soldering. Parts made from the product should not be soldered together with gold or palladium solder. The product is also very suitable for laser welding.

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCMSLQ

FINOFRAME COCR SOLID

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

Safety instructions

Warning: Contains Cobalt (Co). The dust generated during processing of this product may cause cancer if inhaled, may impair fertility and may probably cause genetic defects. It is essential to observe the following safety measures.

Always wear respiratory protection (filter class FFP3), tight-fitting safety goggles, protective gloves and protective clothing and always switch on extraction equipment with filter class Hepa H. Do not inhale dust. Avoid contact with skin, mouth, eyes and clothing. Do not eat or drink while working. Keep away from food and drink. Wash hands after use. Remove contaminated clothing and protective equipment before entering eating areas. Keep away from sources of ignition. Do not smoke.

Storage

No special storage conditions. Store in the original packaging.

Waste disposal

Dispose of the product and packaging in accordance with local/ regional/ national/ international regulations. Do not dispose of with household waste. Do not allow to enter waterways, groundwater or sewage system.

Note

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The Summary of Safety and Clinical Performance according to (EU) MDR 2017 / 745 Article 32 (SSCP) is available in the European database for medical devices (Eudamed) under the associated base UDI-DI as soon as the corresponding database module exists.

<https://ec.europa.eu/tools/eudamed>

Explanation of the markings on the packaging

Confirmation: The product complies with the applicable European directives

Catalog number/ item number

Batch designation


Clear product identification

Number of products in the packaging

Do not reuse

Observe the safety instructions in the operating instructions and instructions for use

www.wp-ifu.de/fino/cocr-em

Path to the electronic instructions for use

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCMSLQ

FINOFRAME COCR SOLID

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

Instructions for use in paper form available within 7 days

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