SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

FINOFRAME ZR / FINOFRAME ZR HT

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

	<b>Gebrauchsanweisung</b> .....	2
	<b>Instructions for use</b> .....	6
	<b>Manuel d'utilisation</b> .....	10
	<b>Istruzioni per l'uso</b> .....	14
	<b>Instrucciones de uso</b> .....	18
	<b>Instrukcje użytkowania</b> .....	22
	<b>Instruktioner för användning</b> .....	26
	<b>Instruktioner til brug</b> .....	30
	<b>Bruksanvisning</b> .....	34
	<b>Käyttöohjeet</b> .....	38
	<b>Gebruiksaanwijzing</b> .....	42
	<b>Οδηγίες χρήσης</b> .....	46
	<b>Upute za uporabu</b> .....	50
	<b>Instrucțiuni de utilizare</b> .....	54
	<b>Naudojimo instrukcijos</b> .....	58
	<b>Návod k použití</b> .....	62
	<b>Návod na použitie</b> .....	66
	<b>Használati utasítás</b> .....	70

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

FINOFRAME ZR / FINOFRAME ZR HT

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

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### Produkt/ Produkttyp

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## Instructions for use

### Product/ Product type

Pre-sintered zirconium dioxide blanks for the fabrication of customised dental restorations

### Material type

ZrO<sub>2</sub>

Class IIa medical device

Ceramic type 2, class 5 according to DIN EN ISO 6872

### Users

Professional users who fabricate customised dental restorations.

### Target group

Patients of all ages and genders who require customised dental restorations

### Indication / Intended use

The product is a medical device intended for the manufacture of dental prostheses intended for long-term use, which are partially inserted into the human body by means of a clinical procedure, but in this case are inserted into the teeth or attached to the tooth substance with suitable fixing material and are therefore to be assigned to class IIa.

Basic frameworks for crowns, bridges, etc. as well as fully anatomical restorations for the dental technique can be milled. The restorations can be coloured with the appropriate FINOFRAME colouring liquids, veneered with ceramic and/or glazed.

### FINOFRAME ZR, FINOFRAME ZR HT:

Crowns and bridges up to 16 units with max. 2 pontics between 2 crowns in the posterior region, with max. 4 pontics between 2 crowns in the anterior region, veneers, inlays, onlays, primary telescopes, bar constructions

### Contraindication

Do not use in case of proven hypersensitivity to one or more ingredients.

Do not use if there is insufficient space.

### Material properties/ Technical data

Values for the material processed according to the instructions for use after final sintering.

Chemical Composition in %	FINOFRAME ZR	
	white	white
ZrO <sub>2</sub> + HfO <sub>2</sub>	94,23 – 94,9	94,23 – 94,9
Y <sub>2</sub> O <sub>3</sub>	4,95 – 5,35	4,95 – 5,35
Al <sub>2</sub> O <sub>3</sub>	0,15 – 0,35	0,15 – 0,35
SiO <sub>2</sub>	0 – 0,02	0 – 0,02
Fe <sub>2</sub> O <sub>3</sub>	0 – 0,01	0 – 0,01
Na <sub>2</sub> O	0 – 0,04	0 – 0,04

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

FINOFRAME ZR / FINOFRAME ZR HT

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

Physical properties	FINOFRAME ZR	
	white	white
Density g/cm <sup>3</sup>	~ 6,05	~ 6,05
Flexural strength MPa	~ 1400	~ 1400
Coefficient of thermal expansion (25°C -500°C) x 10 <sup>-6</sup> /K	10,4 +/- 0,5	10,4 +/- 0,5
Vickers hardness (HV10)	~ 1250	~ 1250
Chemical solubility µg/cm <sup>2</sup>	≤ 100	≤ 100

Chemical Composition in %	FINOFRAME ZR HT			
	-	-	Multi	Dimension
	white	0M1 – 0M3 A1 – D4	A1 – D4	A1 – D4
ZrO <sub>2</sub> + HfO <sub>2</sub>	92,71 - 92,87	92,70 - 93,17	92,70 - 93,17	90,49 - 93,17
Y <sub>2</sub> O <sub>3</sub>	6,93 - 6,97	6,41 - 6,88	6,41 - 6,88	6,41 - 9,21
Al <sub>2</sub> O <sub>3</sub>	0,04 - 0,06	0,04 - 0,06	0,04 - 0,06	0,04 - 0,06
SiO <sub>2</sub>	0 - 0,02	0 - 0,02	0 - 0,02	0 - 0,02
Fe <sub>2</sub> O <sub>3</sub>	0 - 0,01	0,13 - 0,16	0,13 - 0,16	0,01 - 0,16
ER <sub>2</sub> O <sub>3</sub>	---	0 - 0,49	0 - 0,49	0 - 0,78
Co <sub>3</sub> O <sub>4</sub>	---	0 - 0,01	0 - 0,01	0 - 0,01

Physical properties	FINOFRAME ZR HT			
	-	-	Multi	Dimension
	white	0M1 – 0M3 A1 – D4	A1 – D4	A1 – D4
Density g/cm <sup>3</sup>	~ 6,07	~ 6,07	~ 6,07	~ 6,07
Flexural strength MPa	~ 1100	~ 1100	~ 1100	~ 600 (cutting range) to ~ 1100 (neck/dentin area) *
Coefficient of thermal expansion (25°C -500°C) x 10 <sup>-6</sup> /K	10,5 +/- 0,5	10,5 +/- 0,5	10,5 +/- 0,5	10,5 +/- 0,5
Vickers hardness (HV10)	~ 1250	~ 1250	~ 1250	~ 1250
Chemical solubility µg/cm <sup>2</sup>	≤ 100	≤ 100	≤ 100	≤ 100

\*The placement of the restoration within the blank must be observed. See instructions for use.

### Processing

The desired restoration is milled from the selected blank:

The wall thickness for anterior teeth must not be less than 0.5mm and for posterior teeth 0.6mm. The connection point between two splinted anterior crowns must not be less than 8mm<sup>2</sup>. The connection point between an anterior crown and a bridge pontic, two splinted posterior crowns or two anterior bridge pontics must not be less than 9mm<sup>2</sup>. The connection point between a posterior crown and a pontic or two posterior pontics must not be less than 12mm<sup>2</sup>. The cross-section of a bar construction must not be less than 15mm .<sup>2</sup>

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

**FINOFRAME ZR / FINOFRAME ZR HT**

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

With FINOFRAME ZR HT Dimension, the neck/dentin portion always accounts for 50% of the blank height. The connection point must be placed in this area, i.e. always in the lower half of the blank, to ensure sufficient flexural strength.

The dimensions refer to the dimensions after sintering for all product variants.

The enlargement or shrinkage factor is indicated on the blank depending on the milling system used.

The unsintered restorations are finished if necessary. To achieve a perfect surface and maximum translucency, we recommend cleaning the restoration in the FINOHIT ultrasonic unit with distilled water. Please only use our recommended ultrasonic cleaning device. Devices from other companies could be too "strong" and damage the restoration. Swirl the restoration with plastic tweezers in the distilled water of the ultrasonic unit for 5 to 10 seconds until no more "dust clouds" detach from the restoration. Remove excess water from the restoration by blowing it off with oil-free air or by dabbing it dry with cellulose or cotton swabs. Dry the restoration under an infrared lamp or using an oven. The restoration must be completely dry before sintering. Make sure that the water never boils during drying, as this may cause cracks in the zirconium dioxide.

If colouring is desired, this can be done with the corresponding FINOFRAME colouring liquids. Please also refer to the instructions for use. The finished, completely dry restoration is now sintered.

**Sintering**

	<b>Normal programme</b>	<b>Slow Programme</b>	<b>Translucency Programme</b>	<b>Speed Programme</b>
Heating rate	10°C per minute to 900°C	5°C per minute to 900°C	5°C per minute to 900°C	50°C per minute to 1100°C
Holding time	none	20 minutes	none	none
Heating rate	6°C per minute to final temperature	2°C per minute to final temperature	2°C per minute to final temperature	20°C per minute to final temperature
Holding time	at final temperature 90 minutes	at final temperature 120 minutes	at final temperature 120 minutes	at final temperature 30 minutes
Final temperature	1500°C - 1550°C	1500°C - 1550°C	1500°C - 1630°C	1500°C - 1550°C
Cooling down	unregulated in the closed oven	unregulated in the closed oven	unregulated in the closed oven	unregulated in the closed oven

**Reworking**

After sintering, the restoration is adjusted to the working model using the wet grinding technique with diamond-coated grinding instruments, if necessary. Sintered diamonds, corundum stones or tungsten carbide burs should not be used. Overheating must be avoided.

**Firing the ceramic/glaze**

All commercially available firing ceramics/glazes for zirconium dioxide can be used. Please observe the working instructions of the respective manufacturer.

Firing and cooling according to the firing table of the manufacturer of the ceramic material/glaze used.

**Safety instructions**

Warning: Always wear respiratory protection (filter class FFP2), tight-fitting safety goggles, protective gloves and protective clothing and switch on extraction equipment. The dust generated during processing of this product can cause irritation to the skin/eyes/respiratory tract. Avoid inhalation/contact with skin, mouth and eyes. Do not eat or drink while working. Keep away from food and drink. Wash hands after handling. Remove contaminated clothing and protective equipment before entering eating areas. Keep away from sources of ignition. Do not smoke .

**Storage**

Store in a dry place. Protect from moisture. Store in the original packaging.

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

FINOFRAME ZR / FINOFRAME ZR HT

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Fino GmbH • Mangelsfeld 18 • 97708 Bad Bocklet • Germany • +49 9708 909-420

**Waste disposal**

Dispose of the product and packaging in accordance with local/ regional/ national/ international regulations. Do not dispose of with household waste. Do not allow to enter waterways, groundwater or sewage system.

**Note**

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The Summary of Safety and Clinical Performance according to (EU) MDR 2017 / 745 Article 32 (SSCP) is available in the European database for medical devices (Eudamed) under the associated base UDI-DI as soon as the corresponding database module exists.

<https://ec.europa.eu/tools/eudamed>

**Explanation of the labelling on the packaging**

Confirmation: The product complies with the applicable European directives

Catalogue number/ item number

Batch designation

Unique product identification

Number of products in the packaging

Observe the safety instructions in the operating instructions and instructions for use

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Path to the electronic instructions for use

Instructions for use in paper form available within 7 days

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