SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

ARUM SUPREME ML PRO, ARUM HYPERION MD PRO

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Arum 3D Solutions Ltd. • 26C, Lyon Road • Walton-on-Thames KT12 3PU • United Kingdom • +44 20 7388 1390

	Gebrauchsanweisung	2
	Instructions for use	6
	Manuel d'utilisation	10
	Istruzioni per l'uso	14
	Instrucciones de uso	18
	Instrukcje użytkowania	22
	Instruktioner för användning	26
	Instruktioner til brug	30
	Bruksanvisning	34
	Käyttöohjeet	38
	Gebruiksaanwijzing	42
	Οδηγίες χρήσης	46
	Upute za uporabu	50
	Instrucțiuni de utilizare	54
	Naudojimo instrukcijos	58
	Návod k použití	62
	Návod na použitie	66
	Használati utasítás	70

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

ARUM SUPREME ML PRO, ARUM HYPERION MD PRO

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Arum 3D Solutions Ltd. • 26C, Lyon Road • Walton-on-Thames KT12 3PU • United Kingdom • +44 20 7388 1390

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Produkt/ Produkttyp

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Instructions for use

Product/ Product type

Pre-sintered zirconium dioxide blanks for the fabrication of customised dental restorations

Material type

ZrO₂

Class IIa medical device

Ceramic type 2, class 5 according to DIN EN ISO 6872

Users

Professional users who fabricate customised dental restorations.

Target group

Patients of all ages and genders who require customised dental restorations

Indication / Intended use

The product is a medical device intended for the manufacture of long-term dentures, which are partially inserted into the human body through a clinical procedure, but in this case are inserted into the teeth or attached to the tooth substance/abutment with suitable fixing material and are therefore classified as Class IIa.

Both frameworks for crowns, bridges, etc., and fully anatomical restorations for dental technology can be milled. The restorations can be veneered and/or glazed with ceramic.

Crowns and bridges up to 16 units with a maximum of 2 pontics between 2 crowns in the posterior region, with a maximum of 4 pontics between 2 crowns in the anterior region, veneers, inlays, onlays, primary telescopes, bar constructions, implant-supported crowns and bridges, as well as zirconium dioxide superstructures for two-piece abutments or hybrid abutments for dentures.

Contraindication

Do not use in case of proven hypersensitivity to one or more ingredients.

Do not use if there is insufficient space.

Material properties/ Technical data

Values for the material processed according to the instructions for use after final sintering.

ARUM SUPREME ML PRO 0M1 - 0M3 • A1 – D4 Chemical composition in %	
ZrO ₂ + HfO ₂	91,42 – 92,99
Y ₂ O ₃	6,58 – 6,97
Al ₂ O ₃	0,04 – 0,06
SiO ₂	0 – 0,02
Fe ₂ O ₃	0 – 0,23
ER ₂ O ₃	0 – 0,61
Co ₃ O ₄	0 – 0,01
TiO ₂	0 – 0,63

ARUM SUPREME ML PRO 0M1 - 0M3 • A1 – D4 Physical properties	
Density g/cm ³	~ 6,07
Flexural strength MPa	~ 1200
Coefficient of thermal expansion (25°C -500°C)	10,7 +/- 0,5 x 10 ⁻⁶ K ⁻¹
Vickers hardness (HV10)	~ 1250
Chemical solubility µg/cm ²	≤ 100

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

ARUM SUPREME ML PRO, ARUM HYPERION MD PRO

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Arum 3D Solutions Ltd. • 26C, Lyon Road • Walton-on-Thames KT12 3PU • United Kingdom • +44 20 7388 1390

Chemical Composition in %	ARUM HYPERION MD PRO
	0M1 - 0M3 A1 - D4 • Joy 5 - 9
ZrO ₂ + HfO ₂	91,31 – 92,99
Y ₂ O ₃	6,58 – 8,76
Al ₂ O ₃	0 – 0,06
SiO ₂	0 – 0,02
Fe ₂ O ₃	0 – 0,23
ER ₂ O ₃	0 – 0,97
Co ₃ O ₄	0 – 0,01
TiO ₂	0 – 0,63

Physical properties	ARUM HYPERION MD PRO
	0M1 - 0M3 A1 - D4 Joy 5 - 9
Density g/cm ³	~ 6,05
Bending strength MPa	~ 900 (cutting range) to ~ 1200 (Cervical/dentin area)
Coefficient of thermal expansion (25°C - 500°C)	10,3 +/- 0,5 x 10 ⁻⁶ K ⁻¹
Vickers hardness (HV10)	~ 1250
Chemical solubility µg/cm ²	≤ 100

Processing

The desired restoration is milled from the selected blank:

The wall thickness for anterior teeth must not be less than 0.5mm and for posterior teeth 0.6mm. The connection point between two splinted anterior crowns must not be less than 8mm². The connection point between an anterior crown and a bridge pontic, two splinted posterior crowns or two anterior bridge pontics must not be less than 9mm². The connection point between a posterior crown and a pontic or two posterior pontics must not be less than 12mm². The cross-section of a bar construction must not be less than 15mm².

The dimensions refer to the dimensions after sintering for all product variants.

The enlargement or shrinkage factor is indicated on the blank depending on the milling system used.

The unsintered restorations are finished if necessary. In order to achieve a perfect surface and maximum translucency, we recommend cleaning the restoration with distilled water in the "White-Sonic" ultrasonic cleaning device. Please only use our recommended ultrasonic cleaning device. Devices from other companies could be too "strong" and damage the restoration. Swirl the restoration with plastic tweezers in the distilled water of the ultrasonic device for 5 to 10 seconds until no more "dust clouds" detach from the restoration. Remove excess water from the restoration by blowing it off with oil-free air or by dabbing it dry with cellulose or cotton swabs. Dry the restoration under an infrared lamp or using an oven. The restoration must be completely dry before sintering. Make sure that the water never boils during drying, as this may cause cracks in the zirconium dioxide.

The finished, completely dry restoration is now sintered.

Sintering

	Normal Program	Speed Program* (up to 3 links)	Super-Speed Program* (up to 3 links)
Heating rate	10°C/min to 900°C	52°C/min to 1200°C	Super speed sintering under vacuum in the <i>Programat CS6</i> from IVOCLAR. Please follow the manufacturer's instructions.
Holding time	20min	none	
Heating rate	5°C/min to final temperature	6°C/min to 1550°C	
Holding time	120min	30min	
Final temperature	1500°C - 1550°C	1500°C - 1550°C	
Cooling 1	10°C/min to 300°C	55°C/min to 900°C	
Cooling 2	-	-	

*Sintering in the speed process can lead to reduced bending strength.

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

ARUM SUPREME ML PRO, ARUM HYPERION MD PRO

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Arum 3D Solutions Ltd. • 26C, Lyon Road • Walton-on-Thames KT12 3PU • United Kingdom • +44 20 7388 1390

Reworking

After sintering, the restoration is adjusted to the working model using the wet grinding technique with diamond-coated grinding instruments, if necessary. Sintered diamonds, corundum stones or tungsten carbide burs should not be used. Overheating must be avoided.

Firing the ceramic/glaze

All commercially available firing ceramics/glazes for zirconium dioxide can be used. Please observe the working instructions of the respective manufacturer.

Firing and cooling according to the firing table of the manufacturer of the ceramic material/glaze used.

Safety instructions

Warning: Always wear respiratory protection (filter class FFP2), tight-fitting safety goggles, protective gloves and protective clothing and switch on extraction equipment. The dust generated during processing of this product can cause irritation to the skin/eyes/respiratory tract. Avoid inhalation/contact with skin, mouth and eyes. Do not eat or drink while working. Keep away from food and drink. Wash hands after handling. Remove contaminated clothing and protective equipment before entering eating areas. Keep away from sources of ignition. Do not smoke .

Storage

Store in a dry place. Protect from moisture. Store in the original packaging.

Waste disposal

Dispose of the product and packaging in accordance with local/ regional/ national/ international regulations. Do not dispose of with household waste. Do not allow to enter waterways, groundwater or sewage system.

Note

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The Summary of Safety and Clinical Performance according to (EU) MDR 2017 / 745 Article 32 (SSCP) is available in the European database for medical devices (Eudamed) under the corresponding base UDI-DI as soon as the corresponding database module exists.

<https://ec.europa.eu/tools/eudamed>

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

ARUM SUPREME ML PRO, ARUM HYPERION MD PRO

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Arum 3D Solutions Ltd. • 26C, Lyon Road • Walton-on-Thames KT12 3PU • United Kingdom • +44 20 7388 1390

Explanation of the labelling on the packaging

Confirmation: The product complies with the applicable European directives

Catalogue number/ item number

Batch designation

Unique product identification

Number of products in the packaging

Observe the safety instructions in the operating instructions and instructions for use

www.wp-ifu.de/arum/zrpro

Path to the electronic instructions for use

Instructions for use in paper form available within 7 days

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