CE 0483

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCKTLL

CopraTemp

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

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Instructions for use 5

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Bruksanvisning 26

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Instructions for use

Product/ Product type
PMMA blanks for the fabrication of temporary and permanent dental prostheses as well as splints or denture bases.

Material type
Polymethyl methacrylate (PMMA) - Class IIa medical device

Users
Professional users who fabricate customised dental restorations, dental splints or denture bases

Target group
Patients of all ages and genders who require customised dental restorations, dental splints or denture bases

Indication / Intended use
CopraTemp is suitable exclusively for the manufacture of dental products.

CopraTemp and CopraTemp ML
are medical devices intended for the manufacture of dental prostheses for temporary or long-term use, which are partially introduced into the human body through a clinical procedure, in this case inserted into the teeth or attached to the tooth structure with suitable luting material, and are therefore classified as Class IIa.

The tooth-coloured products CopraTemp A1 – D4 and CopraTemp ML 0M1 – 0M3, A1 – D4 are suitable for the fabrication of temporary crowns, bridges, inlays, onlays, etc., which are later replaced by the final dental prosthesis. However, permanent retention in the oral cavity as a permanent restoration is also possible.

CopraTemp Clear, CopraTemp Base and CopraTemp ML Base
are Class IIa medical devices intended for the manufacture of invasive devices, namely occlusal splints, therapeutic splints, bite regulators and surgical guides as well as denture bases, which are intended for long-term use in the oral cavity as far as the pharynx and cannot be absorbed by the mucous membrane.

CopraTemp Clear is colourless and is used to fabricate splints.

The pink-coloured products CopraTemp Base and CopraTemp ML Base are suitable for the fabrication of denture bases.

Permanent retention in the oral cavity is indicated.

Contraindication
Do not use in case of proven hypersensitivity to one or more ingredients.
Do not use if there is insufficient space.

Veneering ceramics/ glaze
All ZrO_2 veneering ceramics/ glazes

Material properties/ Technical data

Material composition
transparent blanks (CopraTemp Clear)

Polymethyl methacrylate	≥ 99%
Residual monomer (MMA)	≤ 1%

pre-coloured blanks (CopraTemp, CopraTemp ML, CopraTemp Base and CopraTemp ML Base)

Polymethyl methacrylate	≥ 98,83%

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Colour pigments	≤ 0,17%
Residual monomer (MMA)	≤ 1%

Mechanical properties of CopraTemp (all product variants)

Flexural strength	> 100 MPa	DIN EN ISO 10477
Flexural strength	> 65 MPa	DIN EN ISO 20795-1/20795-2
Elastic modulus	> 2000 MPa	DIN EN ISO 20795-1/20795-2
Residual monomer	≤ 1%	DIN EN ISO 20795-1/20795-2
Solubility	< 1 µg/mm³	DIN EN ISO 10477
Solubility	< 1,6 µg/mm³	DIN EN ISO 20795-1/20795-2
Water absorption	~ 23 µg/mm³	DIN EN ISO 10477
Water absorption	< 32 µg/mm³	DIN EN ISO 20795-1/20795-2
Vickers hardness	< 20 HV	

Processing
The desired restoration is milled from the selected blank. Milling cutters suitable for PMMA are recommended for processing. The material can be polished using commercially available polishing agents and polishing tools suitable for PMMA. Please ensure that the material is not overheated.

A minimum wall thickness of 1.0 mm must be observed in the occlusal region and 0.6 mm in the cervical region. The connection between crowns or pontics must be at least 10 mm² in the anterior region and at least 15 mm² in the posterior region. No more than 2 pontics are permitted in the posterior region.

Veneering and customisation
CopraTemp and CopraTemp ML can be veneered and customised with suitable resins and composites.
Before applying the veneering material, the surfaces to be veneered must be pretreated in accordance with the instructions of the manufacturer of the veneering system used.
The respective manufacturer's instructions must be observed for surface conditioning, the selection and application of a suitable bonding agent, and the processing and polymerisation of the veneering material.
Only veneering and bonding systems approved by their respective manufacturer for use on PMMA may be used.

Luting of plastic teeth
Plastic teeth can be attached to denture bases made of CopraTemp Base using a suitable luting material. The respective instructions for use of the manufacturer of the bonder, veneering or luting material used must be observed.

Safety instructions
Warning: Always wear respiratory protection (filter class FFP2), tight-fitting safety goggles, protective gloves and protective clothing and switch on extraction equipment. The dust generated during processing of this product can cause irritation to the skin/eyes/respiratory tract. Avoid inhalation/contact with skin, mouth and eyes. Do not eat or drink while working. Keep away from food and drink. Wash hands after handling. Remove contaminated clothing and protective equipment before entering eating areas. Keep away from sources of ignition. Do not smoke.

Storage
Protect from heat and sunlight.

Waste disposal
Dispose of the product and packaging in accordance with local/ regional/ national/ international regulations.


Note
Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The Summary of Safety and Clinical Performance according to (EU) MDR 2017 / 745 Article 32 (SSCP) is available in the European database for medical devices (Eudamed) under the corresponding base UDI-DI as soon as the corresponding database module exists.

https://ec.europa.eu/tools/eudamed

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Explanation of the labelling on the packaging

CE 0483	Confirmation: The product complies with the applicable European directives
REF	Catalogue number/ item number
LOT	Batch designation
UDI	Unique product identification
	Number of products in the packaging
	Please observe the instructions for use
www.wp-ifu.de/PMMA	Path to the electronic instructions for use
	Instructions for use in paper form available within 7 days
MD	Medical device
	Date of manufacture
	Manufacturer
CH REP	Swiss authorised representative
Rx only	Caution: Federal law restricts this product to sale by or on the order of a licensed physician or dentist.

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

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