CE 0483

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCZRMV

CopraSupreme ML PRO, CopraHyperion MD PRO

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

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Instructions for use 6

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Bruksanvisning 34

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Instructions for use

Product/ Product type
Pre-sintered zirconium dioxide blanks for the fabrication of customised dental restorations

Material type
ZrO_2
Class IIa medical device
Ceramic type 2, class 5 according to DIN EN ISO 6872

Users
Professional users who fabricate customised dental restorations.

Target group
Patients of all ages and genders who require customised dental restorations

Indication / Intended use
The product is a medical device intended for the manufacture of dental prostheses for long-term use, which are partially inserted into the human body through a clinical procedure, but are inserted into the teeth or attached to the tooth substance/abudment with suitable luting material and are therefore classified as Class IIa.

Basic frameworks for crowns, bridges, etc., as well as fully anatomical restorations for dental technology can be milled. The restorations can be veneered and/or glazed with ceramic.

Crowns and bridges up to 16 units with max. 2 pontics between 2 crowns in the posterior region, with max. 4 pontics between 2 crowns in the anterior region, veneers, inlays, onlays, primary telescopes, bar constructions, implant-supported crowns and bridges as well as zirconium dioxide abutments for two-part abutments or hybrid abutments for dentures

Contraindication
Do not use in case of proven hypersensitivity to one or more ingredients.
Do not use if there is insufficient space.

Material properties/ Technical data
Values for the material processed according to the instructions for use after final sintering.

CopraSupreme ML PRO 0M1 - 0M3 • A1 – D4 Chemical composition in %	
$ZrO_2 + HfO_2$	91,42 – 92,99
Y_2O_3	6,58 – 6,97
Al_2O_3	0,04 – 0,06
SiO_2	0 – 0,02
Fe_2O_3	0 – 0,23
ER_2O_3	0 – 0,61
Co_3O_4	0 – 0,01
TiO_2	0 – 0,63

CopraSupreme ML PRO 0M1 - 0M3 • A1 – D4 Physical properties	
Density g/cm³	~ 6,07
Flexural strength MPa	~ 1200
Coefficient of thermal expansion (25°C -500°C)	10,7 +/- 0,5 x $10^{-6}K^{-1}$
Vickers hardness (HV10)	~ 1250
Chemical solubility µg/cm²	≤ 100

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Chemical Composition in %	CopraHyperion MD PRO
	0M1 - 0M3 A1 - D4
$ZrO_2 + HfO_2$	91,31 – 92,99
Y_2O_3	6,58 – 8,76
Al_2O_3	0 – 0,06
SiO_2	0 – 0,02
Fe_2O_3	0 – 0,23
ER_2O_3	0 – 0,97
Co_3O_4	0 – 0,01
TiO_2	0 – 0,63

Physical properties	CopraHyperion MD PRO
	0M1 - 0M3 A1 - D4
Density g/cm³	~ 6,05
Bending strength MPa	~ 900 (cutting range) to ~ 1200 (Cervical/dentin area)
Coefficient of thermal expansion (25°C - 500°C)	10,3 +/- 0,5 x 10^{-6} K^{-1}
Vickers hardness (HV10)	~ 1250
Chemical solubility µg/cm²	≤ 100

Processing

The desired restoration is milled from the selected blank:
The wall thickness for anterior teeth must not be less than 0.5mm and for posterior teeth 0.6mm. The connection point between two splinted anterior crowns must not be less than 8mm². The connection point between an anterior crown and a bridge pontic, two splinted posterior crowns or two anterior bridge pontics must not be less than 9mm². The connection point between a posterior crown and a pontic or two posterior pontics must not be less than 12mm². The cross-section of a bar construction must not be less than 15mm².

The dimensions refer to the dimensions after sintering for all product variants.
The enlargement or shrinkage factor is indicated on the blank depending on the milling system used.

The unsintered restorations are finished if necessary. In order to achieve a perfect surface and maximum translucency, we recommend cleaning the restoration with distilled water in the "White-Sonic" ultrasonic cleaning device. Please only use our recommended ultrasonic cleaning device. Devices from other companies could be too "strong" and damage the restoration. Swirl the restoration with plastic tweezers in the distilled water of the ultrasonic device for 5 to 10 seconds until no more "dust clouds" detach from the restoration. Remove excess water from the restoration by blowing it off with oil-free air or by dabbing it dry with cellulose or cotton swabs. Dry the restoration under an infrared lamp or using an oven. The restoration must be completely dry before sintering. Make sure that the water never boils during drying, as this may cause cracks in the zirconium dioxide.

The finished, completely dry restoration is now sintered.

Sintering

	Normal Program	**Speed Program*** (up to 3 links)	**Super-Speed Program*** (up to 3 links)
Heating rate	10°C/min to 900°C	52°C/min to 1200°C	Super speed sintering under vacuum in the *Programat CS6* from IVOCLAR. Please follow the manufacturer's instructions.
Holding time	20min	none	
Heating rate	5°C/min to final temperature	6°C/min to 1550°C	
Holding time	120min	30min	
Final temperature	1500°C - 1550°C	1500°C - 1550°C	
Cooling 1	10°C/min to 300°C	55°C/min to 900°C	
Cooling 2	-	-	

*Sintering in the speed process can lead to reduced bending strength.

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Caution: For a perfect sintering result, we recommend cleaning your sintering furnace regularly, e.g. with our *TA-Clean* cleaning material (observe instructions for use), especially after the previous use of coloring liquids or zirconium dioxide from other manufacturers.

Reworking
After sintering, the restoration is adjusted to the working model using the wet grinding technique with diamond-coated grinding instruments, if necessary. Sintered diamonds, corundum stones or tungsten carbide burs should not be used. Overheating must be avoided.

Firing the ceramic/glaze
All commercially available firing ceramics/glazes for zirconium dioxide can be used. Please observe the working instructions of the respective manufacturer.
Firing and cooling according to the firing table of the manufacturer of the ceramic material/glaze used.

Safety instructions
Warning: Always wear respiratory protection (filter class FFP2), tight-fitting safety goggles, protective gloves and protective clothing and switch on extraction equipment. The dust generated during processing of this product can cause irritation to the skin/eyes/respiratory tract. Avoid inhalation/contact with skin, mouth and eyes. Do not eat or drink while working. Keep away from food and drink. Wash hands after handling. Remove contaminated clothing and protective equipment before entering eating areas. Keep away from sources of ignition. Do not smoke.

Storage
Store in a dry place. Protect from moisture. Store in the original packaging.

Waste disposal
Dispose of the product and packaging in accordance with local/ regional/ national/ international regulations. Do not dispose of with household waste. Do not allow to enter waterways, groundwater or sewage system.

Note
Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The Summary of Safety and Clinical Performance according to (EU) MDR 2017 / 745 Article 32 (SSCP) is available in the European database for medical devices (Eudamed) under the corresponding base UDI-DI as soon as the corresponding database module exists.

https://ec.europa.eu/tools/eudamed

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Explanation of the labelling on the packaging

CE 0483	Confirmation: The product complies with the applicable European directives
REF	Catalogue number/ item number
LOT	Batch designation
UDI	Unique product identification
	Number of products in the packaging
	Observe the safety instructions in the operating instructions and instructions for use
www.wp-ifu.de/zrpro	Path to the electronic instructions for use
7d.	Instructions for use in paper form available within 7 days
MD	Medical device
	Date of manufacture
	Manufacturer
CH REP	Swiss authorised representative
Rx only	Caution: Federal law restricts this product to sale by or on the order of a licensed physician or dentist.

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

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© copyright Whitepeaks Dental Solutions GmbH, page 1 of 73, revision 1